



PIPseq™ UDI-96 Kit

User Guide

Doc ID: FB0003694

Revision: 2.6

State: Release

Table of Contents

| | |
|---|-----------|
| Introduction | 3 |
| Recommended Use | 3 |
| Compatible PIPseq 3' Single Cell RNA Kits | 3 |
| Index Plate Overview | 4 |
| PIPseq Single Cell Gene Expression Library Overview | 5 |
| PIPseq UDI-g6 Kit Protocol | 6 |
| Reagent Preparation | 6 |
| Sample Index PCR (< 96 samples) | 6 |
| Index Pooling Guidelines | 8 |
| Pooling Recommendations | 11 |
| Oligonucleotide Sequences | 13 |
| Document Revision Summary | 14 |
| Legal Notices | 14 |
| Support | 14 |

Introduction

The Fluent BioSciences PIPseq™ UDI-96 Kit enables multiplex sample preparation of PIPseq Single Cell 3' gene expression libraries for next-generation sequencing on the Illumina platform. Each kit provides 96 pre-mixed unique pairs of i5 and i7 index primers packaged in a single-use 96-well breakable plate with a pierceable foil seal.

Recommended Use

The recommended application for this kit is PIPseq library preparation of Single Cell 3' gene expression libraries to be sequenced on an Illumina platform. Users can break away a single column of 8 pre-mixed, unique indexes without thawing the entire plate for each experiment. Further, the pierceable foil seal ensures single use of each pre-mixed index pair to avoid index contamination across wells. The volumes provided are sufficient for the preparation of up to 96 reactions. **The 96-well plate should be stored at -20°C.**

The PIPseq UDI-96 Kit should only be used with Fluent BioSciences PIPseq 3' Single Cell RNA Kits chemistry versions 4.0PLUS and V. These indexes have not been tested for compatibility with other single-cell library preparation kits or technologies.

Note: Users should not mix index primers provided by the PIPseq v4.0PLUS 3' Single Cell RNA Kits and the index primers in the UDI-96 Kit, as the lengths of the sample indexes differ. Instead, it is recommended to utilize index primers from one kit or the other, as determined by the user's desired number of samples for multiplexing.

Users should refer to the kit-specific User Guide before using the PIPseq UDI-96 Kit primers for PIPseq library preparation. Users should verify they have the latest User Guide revision downloaded at www.fluentbio.com.


Compatible PIPseq 3' Single Cell RNA Kits

| Product Name | Catalog Number | User Guide |
|---|----------------------|----------------------|
| PIPseq T2 3' Single Cell RNA Kit v4.0PLUS | FBS-SCR-T2-8-V4.05 | FB0001026 Rev 12.11+ |
| PIPseq T20 3' Single Cell RNA Kit v4.0PLUS | FBS-SCR-T20-4-V4.05 | FB0002130 Rev 8.9+ |
| PIPseq T100 3' Single Cell RNA Kit v4.0PLUS | FBS-SCR-T100-2-V4.05 | FB0003657 Rev 2.9+ |
| PIPseq V T2 3' Single Cell RNA Kit Bundle | FB0005384 | FB0005260 Rev 1.5+ |
| PIPseq V T10 3' Single Cell RNA Kit Bundle | FB0005385 | FB0004762 Rev 1.5+ |
| PIPseq V T20 3' Single Cell RNA Kit Bundle | FB0005386 | FB0005261 Rev 1.5+ |
| PIPseq V T100 3' Single Cell RNA Kit Bundle | FB0005387 | FB0005262 Rev 1.5+ |

Index Plate Overview

The PIPseq UDI-96 plate provides 96 pre-mixed unique pairs of i5 and i7 index primers that are compatible with the PIPseq 3' Single Cell RNA Kits versions 4.0PLUS and V. These pairs are sequentially labeled by column (Figure 1).

NOTE: Index primers provided in the PIPseq V 3' Single Cell RNA Kit are identical to column 3 of the UDI-96 plate. Users should be sure to use indexes from the other columns if using both index sets together for their pooled libraries.

PIPseq Single Cell UDI-96 Kit | FBS-SCR-UDI-96 

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|---|----|----|----|----|----|----|----|----|----|----|----|
| A | 1 | 9 | 17 | 25 | 33 | 41 | 49 | 57 | 65 | 73 | 81 | 89 |
| B | 2 | 10 | 18 | 26 | 34 | 42 | 50 | 58 | 66 | 74 | 82 | 90 |
| C | 3 | 11 | 19 | 27 | 35 | 43 | 51 | 59 | 67 | 75 | 83 | 91 |
| D | 4 | 12 | 20 | 28 | 36 | 44 | 52 | 60 | 68 | 76 | 84 | 92 |
| E | 5 | 13 | 21 | 29 | 37 | 45 | 53 | 61 | 69 | 77 | 85 | 93 |
| F | 6 | 14 | 22 | 30 | 38 | 46 | 54 | 62 | 70 | 78 | 86 | 94 |
| G | 7 | 15 | 23 | 31 | 39 | 47 | 55 | 63 | 71 | 79 | 87 | 95 |
| H | 8 | 16 | 24 | 32 | 40 | 48 | 56 | 64 | 72 | 80 | 88 | 96 |

Figure 1. Plate map for the UDI-96 kit where each well contains a unique i5/i7 index pair and the wells are labeled sequentially by column. This configuration is intended to allow users to break apart the plate by column so that index pairs can be used most efficiently one row at a time.

PIPseq Single Cell Gene Expression Library Overview

The PIPseq 3' Single Cell gene expression library comprises the standard Illumina paired-end sequences P5 and P7. Read 1 contains the barcode information and is distinct between PIPseq chemistry versions while Read 2 contains the gene expression information. **PIPseq utilizes a unique mixture of Nextera and TruSeq sequences for Read 1 and Read 2**, and these libraries are compatible with standard Illumina sequencing. It is recommended to maximize the use of the flow cell with longer reads when possible.

PIPseq version 4.0PLUS: the libraries prepared in this protocol using the UDI-96 Kit for Illumina are **dual-indexed with 10-base i5 and i7 indexes**. **Read 1 length must be ≥ 54 bases** and the recommended **Read 2 length ≥ 63 bases (Figure 2)**.




Figure 2. Structure of the PIPseq single-cell gene expression library for 4.0 and 4.0PLUS libraries prepared using the UDI-96 kit with 10 bp indices.

PIPseq version V: Read 1 length must be ≥ 45 bases and the recommended read 2 length ≥ 72 bases (Figure 3). It is recommended to maximize the use of the flow cell with longer reads when possible.

Please note that the index primers provided in the PIPseq V 3' Single Cell RNA Kit are identical to column 3 of the UDI-96 plate. Users should be sure to use indexes from the other columns if using both index sets together for their pooled libraries.




Figure 3. Structure of the PIPseq single-cell gene expression library for version V chemistry.

PIPseq UDI-96 Kit Protocol

Users should refer to the PIPseq kit-specific User Guide before using the PIPseq UDI-96 Kit primers for PIPseq library preparation. The protocol presented below is an example sample indexing PCR that should replace the corresponding section within the PIPseq kit-specific User Guide.

Reagent Preparation

Upon receipt of the PIPseq UDI-96 Kit, store it in a -20°C freezer.

Before each experiment, determine the number of indexes necessary for sample multiplexing (see Index Pooling Guidelines). To break apart a column (strip), remove the plate from the -20°C freezer, securely grip the desired strip and press along the perforation while pulling the strip away from the rest of the plate. Return the remainder of the plate to the -20°C freezer. Thaw the index strip for 10 min at room temperature before use.

Sample Index PCR (< 96 samples)



Before starting the Sample Index PCR, choose the appropriate sample index sets to ensure that index combinations are color balanced and that no sample index combinations overlap in a multiplexed sequencing run (Index Pooling Guidelines for recommendations). There are 96 pre-mixed i5 and i7 indexes provided with this kit to allow for unique dual indexing of 96 samples.

1. Determine the number of libraries that will be amplified and pooled for sequencing while ensuring that a valid index combination is chosen according to color balance guidelines in the Index Pooling Guidelines section below.
2. Thaw the selected number of index strips at room temperature for 10 minutes.
3. Mix briefly by vortexing and then centrifuge the index strip on a benchtop minifuge to collect all of the primers at the bottom of the well.
4. Thaw 4X PCR Master Mix (alternatively, (HC) Library Prep Mix B) on ice. Once thawed, flick the tubes several times, pipette mix 10 times, and then briefly centrifuge to collect (**DO NOT VORTEX**).
5. Orient the strip according to the provided plate map (Fig 1) to ensure the correct indexes will be removed. With a pipette tip, pierce the desired wells and transfer the volume of primer mix required for the PCR reaction (see step 6).



CAUTION: Ensure pipette tips are not reused across wells to avoid cross-contamination of indexed primers. Do not reuse primer if the seal has been previously pierced to avoid cross-contamination of indexed primers.

6. To the 32.5 μ L of each PIPseq gene expression library (volume depends on PIPseq protocol), add the following reagents in the order in which they appear in the table(s) below. **CAUTION: Please review which kit you have purchased to confirm which table is appropriate.**

PIPseq V Protocol

| Reagent | Volume Per Sample (μ L) |
|--|------------------------------|
| Cleaned library DNA | 32.5 |
| UDI-XX (one well of UDI-96 plate) | 5 |
| 4X PCR Master Mix | 12.5 |
| Total | 50 |

PIPseq V4.0PLUS Protocol

| Reagent | Volume Per Sample (μ L) |
|--|------------------------------|
| Cleaned library DNA | 32.5 |
| UDI-XX (one well of UDI-96 plate) | 5 |
| HC Library Prep Mix B (V4.0PLUS) | 12.5 |
| Total | 50 |

7. Mix the reactions by pipetting up and down 10 times at the 25 μ L stroke, and briefly centrifuge to collect all liquid at the bottom of the tubes.
8. Proceed with the sample index PCR according to the specific PIPseq Kit specific User Guide.

Index Pooling Guidelines

The PIPseq 3' Single Cell gene expression libraries may be pooled for sequencing, taking into account the differences in cell number and read depth requirements. Refer to Illumina documentation as the primary source for appropriate color balance combinations for the selected sequencing platform. **Fluent-validated index combinations are listed in the “Pooling Recommendations” section below.** These libraries are dual-indexed with 10-base i5 and i7 indexes. Sequencing depth will vary based on your application needs but it is recommended to start with a depth of 20,000 reads per input cell.

Users should note that PIPseq index sequences provided with v4.0PLUS 3' Single Cell RNA kits are dual-indexed with 8-base i5 and i7 indexes. This is distinct from the 10-base i5 and i7 indexes provided with PIPseq V 3' Single Cell RNA kits and the PIPseq UDI-96 Kit. Refer to the table below for the read length requirements when using the 10-base indexes provided in the UDI-96 kit. Do not use UDI-17 - UDI-24 from the UDI-96 kit when combining with libraries that were prepared using the indexes provided in the PIPseq V 3' Single Cell RNA kits.

| PIPseq chemistry | Read 1 Length | Read 2 Length | i5/i7 index adapter length |
|------------------|---------------|---------------|----------------------------|
| v4.0PLUS | Min of 54 | Min of 63 | 10 |
| V | Min of 45 | Min of 72 | 10 |

Users may download a sample sheet with the index sequences for use with Illumina Experiment Manager from the PIPseq UDI-96 Kit product page at www.fluentbio.com.

| Well position | Index ID | i5 bases in adapter sequence | i7 bases in adapter sequence |
|---------------|----------|------------------------------|------------------------------|
| A1 | UDI-1 | TCGTGGAGCG | CGCTCAGTTC |
| B1 | UDI-2 | CTACAAGATA | GAATTGAGTG |
| C1 | UDI-3 | TATAGTAGCT | ATATGAGACG |
| D1 | UDI-4 | TGCCTGGTGG | CTTATGGAAT |
| E1 | UDI-5 | ACATTATCCT | TAATCTCGTC |
| F1 | UDI-6 | GTCCACTTGT | GCGCGATGTT |
| G1 | UDI-7 | TGGAACAGTA | AGAGCACTAG |
| H1 | UDI-8 | CCTTGTTAAT | TGCCTTGATC |
| A2 | UDI-9 | GTTGATAGTG | CTACTCAGTC |
| B2 | UDI-10 | ACCAGCGACA | TTCTACAGAA |
| C2 | UDI-11 | CATACACTGT | GAACATACGG |
| D2 | UDI-12 | GTGTGGCGCT | CCTATGACTC |
| E2 | UDI-13 | ATCACGAAGG | TAATGGCAAAG |
| F2 | UDI-14 | CGGCTCTACT | GTGCCGCTTC |

| Well position | Index ID | i5 bases in adapter sequence | i7 bases in adapter sequence |
|---------------|----------|------------------------------|------------------------------|
| G2 | UDI-15 | GAATGCACGA | CGGCAATGGA |
| H2 | UDI-16 | AAGACTATAG | GCCGTAACCG |
| A3 | UDI-17 | TCGGCAGCAA | AACCATTCTC |
| B3 | UDI-18 | CTAATGATGG | TCCAATTCTA |
| C3 | UDI-19 | GGTTGCCTCT | CTAATGATGG |
| D3 | UDI-20 | CGCACATGGC | TCGGCCTATC |
| E3 | UDI-21 | GGCCTGTCTT | AGTCAACCAT |
| F3 | UDI-22 | CTGTGTTAGG | GAGCGCAATA |
| G3 | UDI-23 | TAAGGAACGT | AACAAGGCGT |
| H3 | UDI-24 | CTAACTGTAA | GTATGTAGAA |
| A4 | UDI-25 | GGCGAGATGG | TTCTATGGTT |
| B4 | UDI-26 | AATAGAGCAA | CCTCGCAACC |
| C4 | UDI-27 | TCAATCCATT | TGGATGCTTA |
| D4 | UDI-28 | TCGTATGCGG | ATGTCGTGGT |
| E4 | UDI-29 | TCCGACCTCG | AGAGTGCGGC |
| F4 | UDI-30 | CTTATGGAAT | TGCCTGGTGG |
| G4 | UDI-31 | GCTTACGGAC | TGCGTGTAC |
| H4 | UDI-32 | GAACATACGG | CATACACTGT |
| A5 | UDI-33 | GTCGATTACA | CGTATAATCA |
| B5 | UDI-34 | ACTAGCCGTG | TACGCGGCTG |
| C5 | UDI-35 | AAGTTGGTGA | GCGAGTTACC |
| D5 | UDI-36 | TGGCAATATT | TACGGCCGGT |
| E5 | UDI-37 | GATCACCGCG | GTCGATTACA |
| F5 | UDI-38 | TACCATCCGT | CTGTCTGCAC |
| G5 | UDI-39 | GCTGTAGGAA | CAGCCGATTG |
| H5 | UDI-40 | CGCACTAATG | TGACTACATA |
| A6 | UDI-41 | GACAACCTGAA | ATTGCCGAGT |
| B6 | UDI-42 | AGTGGTCAGG | GCCATTAGAC |
| C6 | UDI-43 | TTCTATGGTT | GGCGAGATGG |
| D6 | UDI-44 | AATCCGGCCA | TGGCTCGCAG |
| E6 | UDI-45 | CCATAAGGTT | TAGAATAACG |
| F6 | UDI-46 | ATCTCTACCA | CAGTAGTTGT |
| G6 | UDI-47 | CGGTGGCGAA | TATCCAGGAC |
| H6 | UDI-48 | TAACAATAGG | AGTGCCACTG |

| Well position | Index ID | i5 bases in adapter sequence | i7 bases in adapter sequence |
|---------------|----------|------------------------------|------------------------------|
| A7 | UDI-49 | CTGGTACACG | GGCCATCATA |
| B7 | UDI-50 | TCAACGTGTA | ACATGGGTGTC |
| C7 | UDI-51 | ACTGTTGTGA | GACAGACAGG |
| D7 | UDI-52 | GTGCGTCCTT | TCTTACATCA |
| E7 | UDI-53 | AGCACATCCT | TTACAATTCC |
| F7 | UDI-54 | TTCCGTCGCA | AAGCTTATGC |
| G7 | UDI-55 | CTTAACCACT | TATTCCTCAG |
| H7 | UDI-56 | GCCTCGGATA | CTCGTGCGTT |
| A8 | UDI-57 | CGTCGACTGG | TTAGGATAGA |
| B8 | UDI-58 | TACTAGTCAA | CCGAAGCGAG |
| C8 | UDI-59 | ATAGACCGTT | GGACCAACAG |
| D8 | UDI-60 | ACAGTTCCAG | TTCCAGGTAA |
| E8 | UDI-61 | AGGCATGTAG | TGATTAGCCA |
| F8 | UDI-62 | GCAAGTCTCA | TAACAGTGTT |
| G8 | UDI-63 | TTGGCTCCGC | ACCGCGCAAT |
| H8 | UDI-64 | AACTGATACT | GTTCGCGCCA |
| A9 | UDI-65 | GTAAGGCATA | AGACACATTA |
| B9 | UDI-66 | AATTGCTGCG | GCGTTGGTAT |
| C9 | UDI-67 | TTACAATTCC | GATAACAAGT |
| D9 | UDI-68 | AACCTAGCAC | TTGTTCCGTG |
| E9 | UDI-69 | TCTGTGTGGA | AAGTACTCCA |
| F9 | UDI-70 | GGAATTCCAA | ACGTCAATAC |
| G9 | UDI-71 | AAGCGCGCTT | GGTGTACAAG |
| H9 | UDI-72 | TGAGCGTTGT | CCACCTGTGT |
| A10 | UDI-73 | ATCATAGGCT | GTTCCGCAGG |
| B10 | UDI-74 | TGTTAGAAGG | ACCTTATGAA |
| C10 | UDI-75 | GATGGATGTA | CGCTGCAGAG |
| D10 | UDI-76 | ACGGCCGTCA | GTAGAGTCAG |
| E10 | UDI-77 | CGTTGCTTAC | GGATAACCAGA |
| F10 | UDI-78 | TGACTACATA | CGCACTAATG |
| G10 | UDI-79 | CGGCCTCGTT | TCCTGACCGT |
| H10 | UDI-80 | CAAGCATCCG | CTGGCTTGCC |
| A11 | UDI-81 | TCGTCTGACT | ACCAGCGACA |
| B11 | UDI-82 | CTCATAGCGA | TTGTAAACGGT |

| Well position | Index ID | i5 bases in adapter sequence | i7 bases in adapter sequence |
|---------------|----------|------------------------------|------------------------------|
| C11 | UDI-83 | AGACACATTA | GTAAGGCATA |
| D11 | UDI-84 | GCGCGATGTT | GTCCACTTGT |
| E11 | UDI-85 | CATGAGTACT | TTAGGTACCA |
| F11 | UDI-86 | ACGTCAATAC | GGAATTCCAA |
| G11 | UDI-87 | GATACCTCCT | CATGTAGAGG |
| H11 | UDI-88 | ATCCGTAAGT | TACACGCTCC |
| A12 | UDI-89 | CGTGTATCTT | GCTTACGGAC |
| B12 | UDI-90 | GAACCATGAA | CGCTTGAAGT |
| C12 | UDI-91 | GGCCATCATA | CGCCTTCTGA |
| D12 | UDI-92 | ACATACTTCC | ATACCAACGC |
| E12 | UDI-93 | TATGTGCAAT | CTGGATATGT |
| F12 | UDI-94 | GATTAAGGTG | CAATCTATGA |
| G12 | UDI-95 | ATGTAGACAA | GGTGAATAC |
| H12 | UDI-96 | CACATCGGTG | TGGACGGAGG |

Pooling Recommendations

Fluent has validated *in silico* the following index combinations in the PIPseq UDI-96 Kit for color balance according to Illumina guidelines for the NextSeq2000 and the NovaSeqX. Not all combinations have been experimentally validated.

| NextSeq 2000 | |
|--------------|--|
| Plexity | Recommended well pairs |
| 2 | A1,B1 : C1,D1 : E1,F1 : G1,H1 (Adjacent well pairs from top) |
| | G2,H2 |
| | C3,D3 |
| | G3,H3 |
| | A4,B4 |
| | E4,F4 |
| | G4,H4 |
| | A5,B5 |
| | C5,D5 |
| | A6,B6 |
| | C6,D6 |
| | E6,F6 |
| | A7,B7 |
| | C7,D7 |

| NextSeq 2000 | |
|-----------------------------------|--|
| Plexity | Recommended well pairs |
| | G7,H7 |
| | A8,B8 |
| | E8,F8 |
| | G8,H8 |
| | A9,B9 |
| | E9,F9 |
| | G9,H9 |
| | A10,B10 |
| | C10,D10 |
| | E10,F10 |
| | A11,B11 |
| | G11,H11 |
| | A12,B12 |
| | G12,H12 |
| 3 | A1,B1,C1 : D1,E1,F1 : F1,G1,H1 (Adjacent wells from top) |
| | COLUMN 2 Adjacent wells from top |
| | COLUMN 3 Adjacent wells from top |
| | COLUMN 4 Adjacent wells from top |
| | A5,B5,C5 |
| | F5,G5,H5 |
| | COLUMN 6 Adjacent wells from top |
| | COLUMN 7 Adjacent wells from top |
| | COLUMN 8 Adjacent wells from top |
| | COLUMN 9 Adjacent wells from top |
| | COLUMN 10 Adjacent wells from top |
| | COLUMN 11 Adjacent wells from top |
| COLUMN 12 Adjacent wells from top | |
| 4 | The first 4 or last 4 wells in any column |
| 5 | The first 5 or last 5 wells in any column |
| 6 | The first 6 or last 6 wells in any column |
| 7 | The first 7 or last 7 wells in any column |
| 8 | The entire column, any column |

| NovaSeq X | |
|-----------|------------------------|
| Plexity | Recommended well pairs |
| 2 | G12,H12 |
| 3 | D2,E2,F2 |
| 3 | D10,E10,F10 |
| 3 | F12,G12,H12 |
| 4 | E1,F1,G1,H1 |

| NovaSeq X | |
|--|--|
| Plexity | Recommended well pairs |
| | A2,B2,C2,D2 |
| | E2,F2,G2,H2 |
| | A3,B3,C3,D3 |
| | E3,F3,G3,H3 |
| | A4,B4,C4,D4 |
| | A5,B5,C5,D5 |
| | E8,F8,G8,H8 |
| | A9,B9,C9,D9 |
| | E9,F9,G9,H9 |
| | A10,B10,C10,D10 |
| | E12,F12,G12,H12 |
| | 5 |
| A4 - E4 | |
| The first 5 or last 5 wells in COLUMN 5 | |
| A6 - E6 | |
| A7 - E7 | |
| The first 5 or last 5 wells in COLUMN 8, 9 or 10 | |
| D12 - H12 | |
| 6 | The first 6 or last 6 wells in COLUMNS 1, 2 or 3 |
| | A4 - F4 |
| | The first 6 or last 6 wells in COLUMN 5 |
| | A6 - F6 |
| 8 | The first 6 or last 6 wells in COLUMNS 7, 8, 9, 10, 11 or 12 |
| 8 | The entire column, any column |

Oligonucleotide Sequences

| Name | Sequence (5' - 3') |
|-----------------|---|
| UDI-## P7 Index | CAAGCAGAAGACGGCATAACGAGATXXXXXXXXXXGTCTCGTGGGCTCGGAGATGT GTATAAGAGACAG |
| UDI-## P5 Index | AATGATACGGCGACCACCGAGATCTACACXXXXXXXXXXACACTCTTTCCTACAG ACGC |

Document Revision Summary

Doc ID: **FB0003694**

Revision: **2.6**

Revision date: **September 2024**

General Changes:

- Updated Pooling recommendations and added PIPseq V compatibility

Legal Notices

© 2024 Fluent BioSciences, Inc (Fluent BioSciences). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of Fluent BioSciences, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product.

Fluent BioSciences may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third party marks or names. The use of products described herein is subject to Fluent BioSciences End User License Agreement, available at www.fluentbio.com/legal-notices, or such other terms that have been agreed to in writing between Fluent BioSciences and the user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

Support

Email: support@fluentbio.com



Fluent BioSciences
150 Coolidge Avenue
Watertown, MA 02472